



## CERTIFICATE OF ANALYSIS

### TRUETROL ADULTERATION CONTROL LEVEL 6

CATALOG NUMBER: A-1326 LOT NUMBER: 071315  
EXP. DATE: 03-2016

PRODUCT NAME: TEST TRUE TRUETROL CONTROL LEVEL 6

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Creatinine 98.0%  
Sodium Nitrite 99.4%  
Pyridinium Chlorochromate 99.4%

#### COLORMETRIC ANALYSIS

USING HITACHI p-Modular:

ANALYTE	UNITS	TARGET	SPECIFICATIONS	RESULTS
pH		11.6	9.86 - 13.34	11.48
Specific Gravity		1.001	1.000 - 1.003	1.0020
Creatinine	(mg/dL)	23	19.55 - 26.45	24.77
Nitrite	(mg/L)	110	93.5 - 126.5	112.30
Chromium (Cr [VI])	(µg/mL)	0	0	0

Kathy A. Chambers

Quality Control

7/13/2015