



CERTIFICATE OF ANALYSIS

TRUETROL ADULTERATION CONTROL LEVEL 6

CATALOG NUMBER: A-1326 LOT NUMBER: 120915
EXP. DATE: 10-2016

PRODUCT NAME: TEST TRUE TRUETROL CONTROL LEVEL 6

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Creatinine 98.0%
Sodium Nitrite 99.4%
Pyridinium Chlorochromate 99.4%

COLORMETRIC ANALYSIS

USING HITACHI p-Modular:

ANALYTE	UNITS	TARGET	SPECIFICATIONS	RESULTS
pH		11.6	11.0 - 12.2	11.43
Specific Gravity		1.001	1.000 - 1.003	1.002
Creatinine	(mg/dL)	23	20.7 - 25.3	24.8
Nitrite	(mg/L)	110	99 - 121	114
Chromium (Cr [VI])	(µg/mL)	0	0	0

Kathy A. Chambers

12/9/2015

Quality Control