



CERTIFICATE OF ANALYSIS

TEST TRUE™ CREATININE CALIBRATOR 200 mg/dL

CATALOG NUMBER:	C-1257 & C-1277	LOT NUMBER:	111016
		EXP. DATE:	11-2017

PRODUCT NAME: TEST TRUE CREATININE CALIBRATOR 200 mg/dL

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Creatinine 98.0%

COLORMETRIC ANALYSIS	TARGET	SPECIFICATIONS	RESULTS
USING HITACHI p-Modular:			
mg/dL	200	180 - 220	199.87

Kathy A. Chambers

 Quality Control

11/10/2016

 Date