



## CERTIFICATE OF ANALYSIS

---

CATALOG NUMBER: O-1351, O-1351L, O-1371/100, O-1371/250      LOT NUMBER: 072617  
EXP. DATE: 08-2018

PRODUCT NAME: Oxidant Reagent, R-1 and R-2

APPEARANCE: R-1 clear, colorless solution      R-2 clear, bluish solution

Total Span (  $\Delta S_2$  -  $\Delta S_1$  )      Specification:  $\geq 1500$       Result: 7972

COLORMETRIC ANALYSIS USING HITACHI MODULAR -P:

QC results satisfactory.

ITEM	UNITS	TARGET	RANGE			RESULTS
Oxidant	mg/L	37.5	33.8	-	41.2	38.1
Oxidant	mg/L	75	67.5	-	82.5	71.6
Oxidant	mg/L	150	135	-	165	150.9
Oxidant	mg/L	200	180	-	220	195
Oxidant	mg/L	300	270	-	330	295

*Kathy A. Chambers*

Quality Control

7/27/2017

Date