



## CERTIFICATE OF ANALYSIS

---

CATALOG NUMBER: 0-1351 LOT NUMBER: 10031114

EXP. DATE: 2015-03

PRODUCT NAME: Oxidant Reagent, R-1 and R-2

APPEARANCE: R-1 clear, colorless solution R-2 clear, bluish solution

Total Span (  $\Delta S_2$  -  $\Delta S_1$  ) Specification:  $\geq 1500$  Result: 2172

COLORMETRIC ANALYSIS USING HITACHI MODULAR -P:

QC results satisfactory.

ITEM	UNITS	TARGET	RANGE	RESULTS
Oxidant	mg/L	37.5	33.8 - 41.2	35.6
Oxidant	mg/L	75	67.5 - 82.5	71.2
Oxidant	mg/L	150	135 - 165	152
Oxidant	mg/L	200	180 - 220	200.2
Oxidant	mg/L	300	270 - 330	334.3

*Charles A. Hay*

Quality Control

3/12/2014

Date