



CERTIFICATE OF ANALYSIS

CATALOG NUMBER: O-1351, O-1351L, O-1371/100, O-1371/250 LOT NUMBER: 100919
EXP. DATE: 10-31-2020

PRODUCT NAME: Oxidant Reagent, R-1 and R-2

APPEARANCE: R-1 clear, colorless solution R-2 clear, bluish solution

Total Span (ΔS_2 - ΔS_1) Specification: ≥ 1500 Result: 6582

COLORMETRIC ANALYSIS USING HITACHI MODULAR -P:

QC results satisfactory.

ITEM	UNITS	TARGET	RANGE			RESULTS
Oxidant	mg/L	37.5	33.8	-	41.2	37.9
Oxidant	mg/L	75	67.5	-	82.5	76.3
Oxidant	mg/L	150	135	-	165	154
Oxidant	mg/L	200	180	-	220	205
Oxidant	mg/L	300	270	-	330	310

Kathy A. Chambers

Quality Control

10/10/2019

Date