



CERTIFICATE OF ANALYSIS

TEST TRUE™ OXIDANT (NITRITE) CALIBRATOR 150 mg/L

CATALOG NUMBER: O-1352 LOT NUMBER: 040115
EXT. EXP. DATE: 10-2016

PRODUCT NAME: TEST TRUE OXIDANT (NITRITE) CALIBRATOR 150 mg/L

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Sodium Nitrite 98.0%

COLORMETRIC ANALYSIS

USING HITACHI p Modular:

UNITS	TARGET	SPECIFICATIONS	RESULTS
mg/L	150	135 - 165	151

Kathy A. Chambers

Quality Control

4/1/2015

Date
