



CERTIFICATE OF ANALYSIS

TEST TRUE™ OXIDANT (NITRITE) CALIBRATOR 150 mg/L

CATALOG NUMBER: O-1352 & O-1372 LOT NUMBER: 061918
EXT. EXP. DATE: 12-2019

PRODUCT NAME: TEST TRUE OXIDANT (NITRITE) CALIBRATOR 150 mg/L

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Sodium Nitrite 98.0%

COLORMETRIC ANALYSIS

USING HITACHI p Modular:

| UNITS | TARGET | SPECIFICATIONS | RESULTS |
|-------|--------|----------------|---------|
| mg/L | 150 | 135 - 165 | 153.0 |

Kathy A. Chambers

Quality Control

6/19/2018

Date
