



CERTIFICATE OF ANALYSIS

TEST TRUE™ OXIDANT (NITRITE) CALIBRATOR 150 mg/L

CATALOG NUMBER: O-1352 & O-1372 LOT NUMBER: 071420
EXT. EXP. DATE: 01-31-2022

PRODUCT NAME: TEST TRUE OXIDANT (NITRITE) CALIBRATOR 150 mg/L

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Sodium Nitrite 98.0%

COLORMETRIC ANALYSIS

USING HITACHI p Modular:

UNITS	TARGET	SPECIFICATIONS	RESULTS
mg/L	150	127.5 - 172.5	149.0

Terra S. Brown

Quality Control

7/22/2020

Date
