



CERTIFICATE OF ANALYSIS

CATALOG NUMBER: P-1058N & P-1078N LOT NUMBER: 092117
EXP. DATE: 10-2018

PRODUCT NAME: TEST TRUE pH 4.0 CALIBRATOR

APPEARANCE: CLEAR, COLORLESS SOLUTION

COLORMETRIC ANALYSIS

USING HITACHI p-Modular:

| | TARGET | SPECIFICATIONS | RESULTS |
|----|--------|----------------|---------|
| pH | 4.0 | 3.8 - 4.2 | 4.00 |

Kathy A. Chambers

Quality Control

9/21/2017

Date
