



CERTIFICATE OF ANALYSIS

SPECIFIC GRAVITY TEST TRUE CONTROL FOR REFRACTOMETER 1.023

CATALOG NUMBER: R-1127 / R-1227 LOT NUMBER: 032614
EXP. DATE: 09-2015

PRODUCT NAME: SPECIFIC GRAVITY CONTROL FOR REFRACTOMETER 1.023

APPEARANCE: CLEAR, COLORLESS SOLUTION

| REFRACTOMETER ANALYSIS | UNITS | TARGET | SPECIFICATIONS | RESULTS |
|---|-------|--------|----------------|---------|
| USING RUDOLPH J57 AUTOMATIC DIGITAL REFRACTOMETER | | | | |
| DIGITAL REFRACTOMETER: | | | | |

| | | | | |
|------------------|--|-------|-------------------|--------|
| Specific Gravity | | 1.023 | 1.01955 - 1.02645 | 1.0230 |
|------------------|--|-------|-------------------|--------|

Charles A. Hay

Quality Control

3/26/2014

Date