



CERTIFICATE OF ANALYSIS

Revised

SPECIFIC GRAVITY TEST TRUE CONTROL FOR REFRACTOMETER 1.023

CATALOG NUMBER: R-1127 / R-1227 LOT NUMBER: 092718
EXP. DATE: 04-2020

PRODUCT NAME: SPECIFIC GRAVITY CONTROL FOR REFRACTOMETER 1.023

APPEARANCE: CLEAR, COLORLESS SOLUTION

REFRACTOMETER ANALYSIS	UNITS	TARGET	SPECIFICATIONS	RESULTS
USING RUDOLPH J57 AUTOMATIC DIGITAL REFRACTOMETER				
DIGITAL REFRACTOMETER:				

Specific Gravity		1.0230	1.0220 - 1.0240	1.0231
------------------	--	--------	-----------------	--------

Kathy A. Chambers

Quality Control

9/27/2018

Date

Revised 10/28/2019