



## CERTIFICATE OF ANALYSIS

---

### TRUETROL ADULTERATION CONTROL LEVEL 1

CATALOG NUMBER: T-3001 & T-3017                      LOT NUMBER: 112519  
 EXP. DATE: 12-31-2020

PRODUCT NAME: TEST TRUE TRUETROL CONTROL LEVEL 1

APPEARANCE: CLEAR, YELLOW SOLUTION

ANALYTE PURITY                      Creatinine                      98.0%  
    Pyridinium Chlorochromate                      99.4%

**COLORMETRIC ANALYSIS**

USING HITACHI p-modular:

	UNITS	TARGET	SPECIFICATIONS	RESULTS
pH	(mg/L)	3.6	3.1 - 4.1	3.6
Specific Gravity	(mg/L)	1.023	1.019 - 1.027	1.023
Creatinine	(mg/dL)	1.25	1.06 - 1.44	1.31
Nitrite	(mg/L)	0	0	0
Chromium	(Cr [VI], µg/mL)	62	53 - 71	61.7

<u>Kathy A. Chambers</u> Quality Control	<u>11/25/2019</u> Date
---	---------------------------

---

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.