



---

## CERTIFICATE OF ANALYSIS

---

### TRUETROL ADULTERATION CONTROL LEVEL 2

CATALOG NUMBER: T-3002 & T-3027                      LOT NUMBER: 041321  
EXP. DATE: 04-30-2022

PRODUCT NAME: TEST TRUE TRUETROL CONTROL LEVEL 2

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY                      Creatinine                      100.0%

#### COLORMETRIC ANALYSIS

USING HITACHI p-Modular:

<u>ITEM</u>	<u>UNITS</u>	<u>TARGET</u>	<u>SPECIFICATIONS</u>	<u>RESULTS</u>
pH		4.25	3.61 - 4.89	4.11
Specific Gravity		1.017	1.014 - 1.020	1.015
Creatinine	(mg/dL)	3.5	2.98 - 4.03	3.70
Nitrite	(mg/L)	0	0	0
Chromium	(Cr [VI], µg/mL)	0	0	0

*Terra S. Brown*

Quality Control

4/19/2021

Date