



CERTIFICATE OF ANALYSIS

TEST TRUE™ CREATININE R1 & R2 REAGENTS

CATALOG NUMBERS:

C-1250 & C-1270

LOT NUMBER: 031325

EXP. DATE: 03-31-2027

PRODUCT NAME:

Creatinine Reagent, R-1 and R-2

APPEARANCE:

R-1 clear, colorless solution

R-2 clear, yellowish solution

COLORMETRIC ANALYSIS:

Using Olympus AU 480

Total Span ($\Delta S2$ - $\Delta S1$) Specification: ≥ 0.038 Result: 0.0499

Linearity Standard Results:

	<u>UNITS</u>	<u>TARGET</u>	<u>SPECIFICATIONS</u>	<u>RESULTS</u>
Creatinine	mg/dL	1.0	0.85 - 1.15	1.01
Creatinine	mg/dL	3.0	2.55 - 3.45	3.04
Creatinine	mg/dL	5.0	4.25 - 5.75	5.05
Creatinine	mg/dL	10.0	8.50 - 11.50	10.06
Creatinine	mg/dL	25.0	21.25 - 28.75	24.70
Creatinine	mg/dL	50.0	42.50 - 57.50	52.13
Creatinine	mg/dL	100.0	85.00 - 115.00	107.00
Creatinine	mg/dL	200.0	170.00 - 230.00	198.14
Creatinine	mg/dL	300.0	255.00 - 345.00	298.11
Creatinine	mg/dL	400.0	340.00 - 460.00	406.02
Creatinine	mg/dL	500.0	425.00 - 575.00	507.39

Control Recovery Results:

	<u>UNITS</u>	<u>TARGET</u>	<u>SPECIFICATIONS</u>	<u>RESULTS</u>
Control 1	mg/dL	1.25	1.06 - 1.44	1.25
Control 2	mg/dL	3.50	2.98 - 4.03	3.38
Control 3	mg/dL	7.00	5.95 - 8.05	7.22
Control 4	mg/dL	23.00	19.55 - 26.45	23.35

Terra S. Brown

Quality Control

3/24/2025

Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.