



CERTIFICATE OF ANALYSIS

TEST TRUE™ CREATININE CALIBRATOR 300 mg/dL

CATALOG NUMBER: C-1258 & C-1278

LOT NUMBER: 040825

EXP. DATE: 04-30-2026

PRODUCT NAME: TEST TRUE CREATININE CALIBRATOR 300 mg/dL

APPEARANCE: CLEAR, COLORLESS SOLUTION

ANALYTE PURITY Creatinine \geq 98.0%

COLORMETRIC ANALYSIS	TARGET	SPECIFICATIONS	RESULTS
USING OLYMPUS AU 480			
mg/dL	300	255 - 345	300.00

Terra S. Brown

Quality Control

4/18/2025

Date

Axiom's calibrators and controls are prepared to target concentrations by gravimetric methods, with the use of balances that are regularly calibrated and verified with traceability to NIST. Critical components in calibrators and controls are traceable to NIST standards. Automated chemistry analyzers used during QA testing are calibrated with or referenced against NIST standards.

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of this product for its particular use. See product insert for additional information.