



CERTIFICATE OF ANALYSIS

TEST TRUE™ NITRITE LINEARITY STANDARD KIT

CATALOG NUMBER: LS5-1010 LOT NUMBER: 112025

EXP. DATE: 2026-05-31

PRODUCT NAME: TEST TRUE NITRITE LINEARITY STANDARD KIT

APPEARANCE: CLEAR, COLORLESS SOLUTIONS

COLORMETRIC ANALYSIS: USING OLYMPUS AU 480

ITEM	TARGET	SPECIFICATIONS	RESULTS
Nitrite	125	106.3 - 143.8	113.25
Nitrite	250	212.5 - 287.5	231.34
Nitrite	500	425 - 575	480.97
Nitrite	1000	850 - 1150	1017.28
Nitrite	1500	1275 - 1725	1585.40

Terra S. Brown
Quality Control

11/20/2025
Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.

