



CERTIFICATE OF ANALYSIS

TEST TRUE™ SPECIFIC GRAVITY LINEARITY STANDARD KIT

CATALOG NUMBER: LS5-1100 & LS6-1110

LOT NUMBER: 091425

EXP. DATE: 2026-03-31

PRODUCT NAME: TEST TRUE SPECIFIC GRAVITY LINEARITY STANDARD KIT

APPEARANCE: CLEAR, COLORLESS SOLUTIONS

INSTRUMENT ANALYSIS: OLYMPUS AU 480 ANALYZER AND RUDOLPH REFRACTOMETER

ITEM	TARGET	SPECIFICATIONS	OLYMPUS RESULTS	REFRACTOMER RESULTS
Specific Gravity	1.001	1.000 - 1.002	1.0008	1.0011
Specific Gravity	1.005	1.004 - 1.006	1.0050	1.0049
Specific Gravity	1.010	1.008 - 1.012	1.0096	1.0102
Specific Gravity	1.020	1.017 - 1.023	1.0200	1.0203
Specific Gravity	1.030	1.025 - 1.035	1.0290	1.0302
Specific Gravity	1.040	1.035 - 1.045	1.0411	1.0424

Terra S. Brown

Quality Control

9/15/2025

Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.