



CERTIFICATE OF ANALYSIS

TEST TRUE™ SPECIFIC GRAVITY LINEARITY STANDARD KIT

CATALOG NUMBER: LS5-1100 & LS6-1110

LOT NUMBER: 101823

EXP. DATE: 04-30-2024

PRODUCT NAME: TEST TRUE SPECIFIC GRAVITY LINEARITY STANDARD KIT

APPEARANCE: CLEAR, COLORLESS SOLUTIONS

INSTRUMENT ANALYSIS: OLYMPUS AU 480 ANALYZER AND RUDOLPH REFRACTOMETER

ITEM	TARGET	SPECIFICATIONS	OLYMPUS RESULTS	REFRACTOMER RESULTS
Specific Gravity	1.001	1.000 - 1.002	1.001	1.0010
Specific Gravity	1.005	1.004 - 1.006	1.005	1.0049
Specific Gravity	1.010	1.008 - 1.012	1.009	1.0099
Specific Gravity	1.020	1.017 - 1.023	1.020	1.0203
Specific Gravity	1.030	1.025 - 1.035	1.028	1.0302
Specific Gravity	1.040	1.035 - 1.045	1.038	1.0429

Terra S. Brown

Quality Control

10/31/2023

Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.