



CERTIFICATE OF ANALYSIS

TEST TRUE™ CHLORIDE SPECIFIC GRAVITY LINEARITY STANDARD KIT

CATALOG NUMBER: LS5-1101CHL LOT NUMBER: 012025
EXP. DATE: 07-31-2025

PRODUCT NAME: TEST TRUE CHLORIDE SPECIFIC GRAVITY LINEARITY STANDARD KIT

APPEARANCE: CLEAR, COLORLESS SOLUTIONS

INSTRUMENT ANALYSIS: OLYMPUS AU 480 ANALYZER

ITEM	TARGET	SPECIFICATIONS	OLYMPUS RESULTS
Specific Gravity	1.002	1.000 - 1.003	1.0015
Specific Gravity	1.005	1.003 - 1.007	1.0051
Specific Gravity	1.010	1.007 - 1.013	1.0097
Specific Gravity	1.020	1.017 - 1.023	1.0204
Specific Gravity	1.030	1.025 - 1.035	1.0309

Terra S. Brown

Quality Control

1/31/2025

Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.