



CERTIFICATE OF ANALYSIS

TEST TRUE™ CHLORIDE SPECIFIC GRAVITY LINEARITY STANDARD KIT

CATALOG NUMBER: LS5-1101CHL LOT NUMBER: 112125
EXP. DATE: 2026-05-31

PRODUCT NAME: TEST TRUE CHLORIDE SPECIFIC GRAVITY LINEARITY STANDARD KIT

APPEARANCE: CLEAR, COLORLESS SOLUTIONS

INSTRUMENT ANALYSIS: OLYMPUS AU 480 ANALYZER

ITEM	TARGET	SPECIFICATIONS	OLYMPUS RESULTS
Specific Gravity	1.002	1.000 - 1.003	1.0028
Specific Gravity	1.005	1.003 - 1.007	1.0060
Specific Gravity	1.010	1.007 - 1.013	1.0101
Specific Gravity	1.020	1.017 - 1.023	1.0197
Specific Gravity	1.030	1.025 - 1.035	1.0292

Terra S. Brown

Quality Control

12/4/2025

Date

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.