



## CERTIFICATE OF ANALYSIS

---

### TEST TRUE™ CHROMATE LINEARITY STANDARD KITS

CATALOG NUMBER:                    LS5-1150 & LS6-1150                    LOT NUMBER: 090923  
EXP. DATE: 03-31-2024

PRODUCT NAME:            TEST TRUE CHROMATE LINEARITY STANDARD KIT

APPEARANCE:                    CLEAR, YELLOW SOLUTIONS

COLORMETRIC ANALYSIS :            USING OLYMPUS AU 480

<u>ITEM</u>	<u>TARGET</u>	<u>SPECIFICATIONS</u>	<u>RESULTS</u>
CHROMATE	10 mg/L	9.0 - 11.0	10.1
CHROMATE	25 mg/L	22.5 - 27.5	25.1
CHROMATE	50 mg/L	45 - 55	48.8
CHROMATE	100 mg/L	90 - 110	98.5
CHROMATE	150 mg/L	135 - 165	146.2
CHROMATE	200 mg/L	180 - 220	193.8

*Terra S. Brown*  
Quality Control

9/29/2023  
Date

---

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.