



## CERTIFICATE OF ANALYSIS

---

### TEST TRUE™ CREATININE LINEARITY STANDARD KIT

CATALOG NUMBERS: LS5-1250, LS5-1260, LS5-1270 & LS5-1300      LOT NUMBER: 071925

EXP. DATE: 2026-01-31

PRODUCT NAME:              TEST TRUE CREATININE LINEARITY STANDARD KIT

APPEARANCE:                      CLEAR, COLORLESS SOLUTIONS

COLORMETRIC ANALYSIS :              USING OLYMPUS AU 480

| <u>ITEM</u> | <u>TARGET</u> | <u>SPECIFICATIONS</u> | <u>RESULTS</u> |
|-------------|---------------|-----------------------|----------------|
| CREATININE  | 1 mg/dL       | .90 - 1.10            | 1.01           |
| CREATININE  | 3 mg/dL       | 2.70 - 3.30           | 3.01           |
| CREATININE  | 5 mg/dL       | 4.50 - 5.50           | 5.01           |
| CREATININE  | 10 mg/dL      | 9.00 - 11.00          | 10.07          |
| CREATININE  | 25 mg/dL      | 22.50 - 27.50         | 24.90          |
| CREATININE  | 50 mg/dL      | 45 - 55               | 53             |
| CREATININE  | 100 mg/dL     | 90 - 110              | 104            |
| CREATININE  | 200 mg/dL     | 180 - 220             | 200            |
| CREATININE  | 300 mg/dL     | 270 - 330             | 300            |
| CREATININE  | 400 mg/dL     | 360 - 440             | 413            |
| CREATININE  | 500 mg/dL     | 450 - 550             | 510            |

*Terra S. Brown*  
Quality Control

7/22/2025  
Date

---

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.