



## CERTIFICATE OF ANALYSIS

---

CATALOG NUMBER: M-1150L, M-1170/500, M-1170/100

LOT NUMBER: 021622

EXP. DATE: 02-29-2024

Reagent Date of Manufacture 2/16/2022

PRODUCT NAME: CHROMATE REAGENT R1

APPEARANCE: CLEAR, COLORLESS PARTICULATE FREE SOLUTION

Total Span ( S2 - S1 ) :      Specification:       $\geq 400$       Result: 717

### COLORMETRIC ANALYSIS

USING HITACHI Modular-p:

<u>ITEM</u>	<u>UNITS</u>	<u>TARGET</u>	<u>SPECIFICATIONS</u>	<u>RESULTS</u>
Chromate	mg/L	10	9.0 - 11.0	<u>10</u>
Chromate	mg/L	25	22.5 - 27.5	<u>25</u>
Chromate	mg/L	50	45 - 55	<u>49</u>
Chromate	mg/L	100	90 - 110	<u>101</u>
Chromate	mg/L	150	135 - 165	<u>150</u>

Terra S. Brown

Quality Control

2/22/2022

Date

---

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional information.