



CERTIFICATE OF ANALYSIS

TEST TRUE™ NITRITE CALIBRATOR 1,000 mg/L

CATALOG NUMBER: N-1012 & N-1072

LOT NUMBER: 101425

EXP. DATE: 2027-04-30

PRODUCT NAME: TEST TRUE NITRITE CALIBRATOR 1,000 mg/L

APPEARANCE: CLEAR, COLORLESS SOLUTION

COLORMETRIC ANALYSIS

USING OLYMPUS AU 480

TARGET	SPECIFICATIONS	RESULTS
1,000	850 -1,150	1,067

Terra S. Brown

Quality Control

10/21/2025

Date

Axiom's calibrators and controls are prepared to target concentrations by gravimetric methods, with the use of balances that are regularly calibrated and verified with traceability to NIST. Critical components in calibrators and controls are traceable to NIST standards. Automated chemistry analyzers used during QA testing are calibrated with or referenced against NIST standards.

Axiom warrants that its products conform to the information contained in this and other Axiom publications. Purchaser must determine the suitability of the product for its particular use. See product insert for additional